

Gojira Fine Chemicals, LLC

A World of Life Science Products - Just A Dragon's Breath Away 9088 Superior Ave. - Suite 202A - Streetsboro, OH 44241

www.gojirafc.com • sales@gojirafc.com

Certificate of Analysis

Pristane

Product #: PR1001

Lot #: 020221

CAS #: 1921-70-6

Synonym: 2,6,10,14-Tetramethylpentadecane

Molecular Formula: $C_{19}H_{40}$

Molecular Weight: 268.5

MFG. Date: Jan. 2012

Test	Specification	Results
Appearance	Pale yellow to clear colorless liquid	Pale yellow to clear colorless liquid
Biological source	Synthetic	Synthetic
Solubility	Soluble in Dichloromethane,	Complies
	Hexane, Toluene, Tetrahydrofuran	
Density at 25°C	0.76 – 0.79 g/ml	0.778 g/ml
Identification by ¹ H-NMR	Should conform to the structure	Conforms to the structure
Chromatographic purity		
by GC (by area	≥ 98%	98.56%
normalization)		

Approved By:

Robert J. Magada
Gojira Fine Chemicals, LLC
Vice President & Technical Director

Tel: 330-974-9505 Fax: 888-211-5523

rob@gojirafc.com